

# EarthResin



## Epoxy Worksurfaces

**FOR SUSTAINABLE DESIGN**




**KEWAUNEE®**

*...encouraging new discovery*


Sustainable  
Epoxy Resin  
Worksurfaces  
and Sinks  
for Labs of the  
21<sup>st</sup> Century

# EarthResin by Kewaunee

## EarthResin Colors



Kewaunee EarthResin has the same physical properties and chemical resistance as Kewaunee Kemresin®.

Printed in the USA on 100% recycled post consumer fiber 

## Kewaunee EarthResin Worksurfaces and Sinks are:

- a full thickness uniform mixture of modified epoxy resin, with a minimum of 10% post consumer recycled content (by weight)
- designed and engineered to meet SEFA performance specifications, assuring laboratories that will last for decades, matching the building life cycle
- manufactured using modern, environmentally responsible methods and materials that meet or exceed all local, state, and national regulations
- manufactured in the United States of America



## Kewaunee Scientific Corporation

- is a member of the U.S. Green Building Council
- is the recipient of the Environmental Excellence Award



## Environmental Statement

Kewaunee Scientific Corporation recognizes that our environment is the responsibility of not only the individual, but the corporate community as well. The success of our customers, our employees, and our company depends on the ability of us all to sustain the resources that underlie the products and services we offer and use. Kewaunee has a two-fold environmental commitment:

1. To provide laboratory furniture, fume hoods, and worksurfaces that minimize energy consumption and adverse environmental impact by designing, engineering, and manufacturing products
  - that use renewable and recyclable resources
  - that use less energy and resources to install and operate, and
  - that match the building life cycle, withstanding the rigors of decades of use, relocation and reuse
2. To manufacture these products in an environmentally responsible manner.



## Kewaunee and LEED Certification

Kewaunee's commitment to the environment includes support of the U.S. Green Building Council's LEED program (Leadership in Energy and Environmental Design). The LEED program is a certification system designed to promote buildings that are environmentally responsible, cost effective, and healthy places for people to live and work.



 **KEWAUNEE®**  
Scientific Corporation

P.O. Box 1842 • Statesville, NC 28687-1842  
Phone: (704) 873-7202 • Fax: (800) 932-3296  
E-mail: [kscmarketing@kewaunee.com](mailto:kscmarketing@kewaunee.com)  
Web Site: [www.kewaunee.com](http://www.kewaunee.com)