

# Kewaunee Steel Cabinet Construction



Kewaunee's fully-welded, heavy-gauge, steel base cabinets are built to rigorous benchmarks. They are certified to meet all SEFA 8 performance standards. Base cabinets are available in many styles of cupboard units, drawer units and drawer and cupboard units. They are available in standard lengths of 12", 15", 18", 21", 24", 30", 36", 42", and 48". Standard heights are 30", 33", and 35-3/4" and standard depth is 22".

## CABINET BODY

- Heavy-gauge formed and welded steel cabinet body construction provides strength and durability.
- Full height removable back (supplied for all cupboard units, and for other units when specified) provides full access to service chase.
- Shelf adjustment and/or drawer channel mounting holes on 1/2" centers.
- Flush corner construction for smooth seamless cabinet face.
- Lipped end panels and rails for built-in dust protection

## STEEL DOOR

- Two-piece door construction painted before assembly. Structurally sound deadened for added rigidity and quiet operation.
- 5-knuckle institutional type stainless steel hinge, 2-1/2" high.
- Continuous hinge reinforcement.
- Positive latching mechanism.
- Rubber door bumpers.
- Variety of door and drawer pulls available.

## STEEL DRAWER

- Full access drawer with lock open feature. Self-closing for safety and housekeeping.
- Interchangeable, interlocking drawer head, sound deadened for quiet operation.
- One piece drawer body with radiused bottom for easy cleaning.
- Radiused drawer suspension provides 150 pound load capacity.
- Nylon tired, stainless steel ball bearing drawer rollers for smooth and quiet operation.

## CABINET BOTTOM & BASE

- One-piece case bottom and horizontal front rail. Adjusting bolt access hole with removable plug button.
- Completely enclosed toe space promotes laboratory cleanliness.
- Lower corner reinforcing gusset with 2-1/2" high leveling bolt (all four corners).

## STEEL FINISH

- Environmentally friendly, electrostatically applied powder coat finish.
- Laboratory grade formulation meets all SEFA 8 chemical requirements.

## ADJUSTABLE SHELF

- Full length and width adjustable shelf for maximum storage.
- Heavy-duty plated steel shelf clips.

## WOOD FINISH

- Synthetic varnish system. Formulation meets all SEFA 8 requirements. - Style 16

## WOOD DOOR AND DRAWER FRONTS

- 3/4" Oak or Maple banded door with hardwood veneer. - Style 16



Exclusive positive door latching mechanism is self-closing and prevents rebounding.



# Kewaunee Steel Cabinet Styles



## Style -01

Square-edge design with inset doors and drawers. Available with a wide selection of pull options. Style 01 is available in both painted and stainless steel.



## Style -02

Square-edge design with inset doors and drawers. Furnished with an anodized aluminum recessed pull.



## Style -11

Square-edge design with full overlay doors and drawers. Available with a wide selection of pull options.



## Style -12

Square-edge design with full overlay doors and drawers. Furnished with an anodized aluminum recessed pull.



## Style -14

Radiused-edge, full overlay doors and drawers. Available with a wide selection of pull options.



## Style -15

Radiused-edge, full overlay doors and drawers. Furnished with an anodized aluminum recessed pull.



## Style -16

Steel cabinet body with wood doors and drawers. Full overlay doors and drawers available with a wide selection of pull options.

