

Wood Cabinet Information

Catalog Number Explanation

E 4 0 W 3 6 2 2 3 6 - 1 0 1 1 AB1C...

<p>Cabinet Type</p> <p>A = Apron Assemblies and Legs B = Open Cupboard Base Cabinet C = All Cupboard Base Cabinet D = All Drawer Base Cabinet E = Drawer & Cupboard Base Cabinet G = Special Purpose Base Cabinet S = Full Height Storage Cabinet W = Wall Storage Cabinet</p> <p>Door & Drawer Configuration</p> <p>Unit Construction and Design</p> <p>Cabinet Height</p> <p>Base Cabinets 31 = 30" High 34 = 32¹/₂" High 36 = 35" High</p> <p>Aprons 3 = 3" High 5 = 4¹⁷/₃₂" High Full Height Cabinets 84 = 84" High</p> <p>Wall Cabinet Height 24 = 24" High 30 = 30" High 36 = 36" High 48 = 48" High</p>	<p>Cabinet Options (See Cabinet Option list below)</p> <p>Style – Hardware – Material Selections (See Selection list below)</p> <p>Special Configurations L = Reverse Construction/Left Hand Swing A = Adjustable Length (Apron Assemblies)</p> <p>Cabinet Length Standard Lengths: 15" – 18" – 21" – 24" 30" – 36" – 42" – 48"</p> <p>Cabinet Depth (Nominal)</p> <p>Base Cabinets 16 = 16" 22 = 22" 28 = 27¹/₂" Aprons 22 = 22" 28 = 28" Wall Cabinets & Full Height Cabinets 12 = 12³/₄" 16 = 16³/₄" 22 = 22"</p>
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Style – Hardware – Material Selections

E 4 0 W 3 6 2 2 3 6 - 1 0 1 1 AB1C...

<p>Door and Drawer Head Style</p> <p>1 = Conventional Radius Lip 3 = Silhouette 5 = Contemporary Full Overlay</p> <p>Hinge Type</p> <p>0 = Stainless Steel – Five Knuckle 2 = Black Powder Coated – Five Knuckle 3 = White Powder Coated – Five Knuckle</p>	<p>Pull Type</p> <p>0 = No Pull * 1 = Rounded Aluminum 2 = Black Wire 3 = White Wire 4 = Stainless Steel Wire 5 = Integral Machined 6 = Chrome Plated Wire 8 = Square Aluminum * pull holes bored on 4" centers but no pull provided</p> <p>Cabinet Exterior</p> <p>1 = Plain Sliced Red Oak 2 = Plain Sliced White Maple †</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

† Maple exterior available for Style 1 and Style 5 only


Cabinet Options – (Add to End of Catalog Number)

- | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>A – Unfinished Cabinet End Panels</p> <p>B1 – 4" Dia. Casters (4) w/ brakes</p> <p>B2 – 4" Dia. Casters (4) w/o brakes</p> <p>B3 – 4" Dia. Casters (2) w/ brakes - (2) w/o brakes</p> <p>B4 – 3" Dia. Casters (4) w/ brakes</p> <p>B5 – 3" Dia. Casters (4) w/o brakes</p> <p>B6 – 3" Dia. Casters (2) w/ brakes - (2) w/o brakes</p> <p>C – 1" Added to Cabinet Height</p> <p>D – Full Top Frame</p> <p>E – Finished Back</p> <p>F – Removable Back (Base cabinet cupboards only)</p> <p>G – Toespace Grill</p> <p>H – Label Holders *</p> <p>J – Disk Locks or J2 – Pin Locks w/Lock Panels*</p> <p>K – Disk Locks or K2 – Pin Locks w/o Lock Panels*</p> | <p>L – Padlock Hasps *</p> <p>M – Levelers for base and full height cabinets</p> <p>O – Shelf Lips</p> <p>P1 – One Pullout Tray - at bottom of cupboard</p> <p>P2 – Two Pullout Trays - center & bottom</p> <p>P5 – One Pullout Pan - at bottom of cupboard</p> <p>P6 – Two Pullout Pans - center & bottom</p> <p>R – 1/8" Tempered Safety Glass</p> <p>T – 100# Full Extension Lineal Ball Bearing Slide</p> <p>W – White Coated Drawer Bottom</p> <p>Y – Intermediate Rails on Style 5 Cabinets</p> <p>U1 – Toe Space on Left End of Base Cabinet</p> <p>U2 – Toe Space on Right End of Base Cabinet</p> <p>U3 – Toe Space on Both Ends of Base Cabinet</p> |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

* Furnished on all cupboards and drawers on cabinet. — Lock keying and/or numbering must be specified.

Wood Cabinet Information

Base Cabinet Dimensional Information



Base cabinets are designed on a 32mm high "unit" system to allow maximum flexibility in the selection of drawer and/or cupboard door sizes.

Standing height cabinets (35" high) have an opening 24 "units" high. Sitting height cabinets (30" high) have an opening 20 "units" high. ADA height cabinets (32 1/2" high) have an opening 22 "units" high. What this means is that any combination of drawer and cupboard door heights equaling 24, 22, or 20 units can be used in a base cabinet.

The adjacent tables indicate the inside clearance dimensions of the base cabinet drawers and cupboards shown in this catalog. Simply match the number code shown on the cabinet to the appropriate code in tables

The dimensions shown in these tables are minimums. They may be slightly larger on some cabinet types and/or styles .

These tables are intended to show drawer and cupboard clearances, not actual construction dimensions. If actual dimensions are required, contact your Kewaunee representative.

Wood Drawer Size Codes

Units High	Drawer Opening Height	Inside Dimensions	
		Depth/Height	Length (Front-to-Back)
2		Pullboard 14"	
3	2 ¹³ / ₁₆ "	2 ¹ / ₈ " †	16 ³ / ₄ "
4	4 ¹ / ₁₆ "	3 ³ / ₈ " †	16 ³ / ₄ "
5	5 ⁵ / ₁₆ "	4 ⁵ / ₈ " †	16 ³ / ₄ "
6	6 ⁹ / ₁₆ "	5 ⁷ / ₈ " †	16 ³ / ₄ " *
7	7 ¹³ / ₁₆ "	7 ¹ / ₈ " †	16 ³ / ₄ "
8	9 ¹ / ₈ "	8 ³ / ₈ " †	16 ³ / ₄ "
10	11 ⁵ / ₈ "	10 ⁷ / ₈ " †	16 ³ / ₄ "
12	14 ¹ / ₈ "	13 ¹ / ₈ " †	16 ³ / ₄ "

† 1/4" less if drawer is top drawer in cabinet

* 11³/₄" on 16" deep base cabinets

Wood Drawer Width I.D. (side-to-side)

Base Cabinet Length	Number of Drawers Across				
	1	2	3	4	6
15"	11 1/2"				
18"	14 1/2"				
21"	17 1/2"				
24"	20 1/2"	8 ³ / ₄ "			
30"	26 1/2"	11 ³ / ₄ "			
36"	32 1/2"	14 ³ / ₄ "	8 ³ / ₄ "		
42"	38 1/2"	17 ³ / ₄ "	10 ³ / ₄ "		
48"	44 1/2"	20 ³ / ₄ "	12 ³ / ₄ "	8 ⁷ / ₈ "	
60"		26 ³ / ₄ "	16 ³ / ₄ "	10 ⁷ / ₈ "	
72"		32 ³ / ₄ "	20 ³ / ₄ "	14 ⁷ / ₈ "	8 ⁷ / ₈ "

Wood Cupboard Size Codes

Units High	Opening Height	Inside Dimensions	
		Height (Top-to-Bottom)	Depth (Front-to-Back)
12	14 ³ / ₈ "	14 ³ / ₈ "	19 ¹ / ₂ "
15	18 ¹ / ₈ "	18 ¹ / ₈ "	19 ¹ / ₂ "
16	19 ³ / ₈ "	19 ³ / ₈ "	19 ¹ / ₂ "
17	20 ⁵ / ₈ "	20 ⁵ / ₈ "	19 ¹ / ₂ "
18	21 ⁷ / ₈ "	21 ⁷ / ₈ "	19 ¹ / ₂ " *
20	24 ¹ / ₄ "	24 ¹ / ₄ "	19 ¹ / ₂ "
22	26 ³ / ₄ "	26 ³ / ₄ "	19 ¹ / ₂ "
24	29 ¹ / ₄ "	29 ¹ / ₄ "	19 ¹ / ₂ "

* 14¹/₂" on 16" deep base cabinets
26" on 27¹/₂" deep base cabinets

Wood Cupboard Width I.D. (side-to-side)

Base Cabinet Length	Number of Cupboards Across				
	1	2	3	4	6
15"	13 1/2"				
18"	16 1/2"				
21"	19 1/2"				
24"	22 1/2"	10 ³ / ₄ "			
30"	28 1/2"	13 ³ / ₄ "			
36"	34 1/2"	16 ³ / ₄ "	10 ³ / ₄ "		
42"	40 1/2"	19 ³ / ₄ "	12 ³ / ₄ "		
48"	46 1/2"	22 ³ / ₄ "	14 ³ / ₄ "	10 ⁷ / ₈ "	
60"		28 ³ / ₄ "	18 ³ / ₄ "	13 ⁷ / ₈ "	
72"		34 ³ / ₄ "	22 ³ / ₄ "	16 ⁷ / ₈ "	10 ⁷ / ₈ "

The dimensions shown in these tables are minimums. They may be slightly larger on some cabinet types and/or styles .

Wood cabinets are available with finished or unfinished end panels
FOR UNFINISHED END PANELS SELECT OPTION A

Unfinished End Panels

When a cabinet's end panel will not be exposed to view the cabinet may be furnished with unfinished ends. The A Option provides unfinished ends on both ends of the cabinet.

Unfinished ends may be specified for base cabinets, full height cabinets or wall cabinets.

Note: Unfinished ends can not be specified on Mobile and Suspended Cabinets, on Center Bench End Sink Cabinets or on cabinets with Toe Spaces on the end of the cabinet.

To specify Unfinished Ends, add **A** to the end of the cabinet number.

Example: **E40D362236-1011A**

Option A