

Wood Utility Tables


Utility Tables are designed for two students and are available in both 31" high sitting and 36" high standing heights. They are also available in 21" and 24" depths, with your choice of a plain apron rail or book/tote tray storage compartments. Tables are available with a 1" thick work top in optional materials.

Adjustable plastic swivel glides are provided on each leg to level table and for easy relocation.

Note: Work tops are available in three materials: Kemresin, Plastic Laminate and Acid Resistant Plastic Laminate. To specify top material, replace the box in part number with the top material designation letter.

Example: T11W312154P0010 specifies an Oak Plain Apron Rail Table with a Plastic Laminate work top.

Utility Tables



Plain Apron Rail Table 20⁷/₈" Deep x 54" Long

Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate
Legs: Adjustable Plastic Swivel Glides

| Stretchers | Height | with Work Top | without Work Top |
|------------|--------|-----------------|------------------|
| No | 31" | T11W312154□00_0 | T11W312154-T00_0 |
| No | 36" | T11W362154□00_0 | T11W362154-T00_0 |
| Yes | 31" | T10W312154□00_0 | T10W312154-T00_0 |
| Yes | 36" | T10W362154□00_0 | T10W362154-T00_0 |





Table with Book Compartments 20⁷/₈" Deep x 54" Long

Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate
Legs: Adjustable Plastic Swivel Glides

| Stretchers | Height | with Work Top | without Work Top |
|------------|--------|-----------------|------------------|
| No | 31" | T12W312154□00_0 | T12W312154-T00_0 |
| No | 36" | T12W362154□00_0 | T12W362154-T00_0 |
| Yes | 31" | T02W312154□00_0 | T02W312154-T00_0 |
| Yes | 36" | T02W362154□00_0 | T02W362154-T00_0 |



Plain Apron Rail Table 24" Deep x 54" Long

Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate
Legs: Adjustable Plastic Swivel Glides

| Stretchers | Height | with Work Top | without Work Top |
|------------|--------|-----------------|------------------|
| No | 31" | T18W312454□00_0 | T18W312454-T00_0 |
| No | 36" | T18W362454□00_0 | T18W362454-T00_0 |
| Yes | 31" | T08W312454□00_0 | T08W312454-T00_0 |
| Yes | 36" | T08W362454□00_0 | T08W362454-T00_0 |




Table with Book Compartments 24" Deep x 54" Long

Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate
Legs: Adjustable Plastic Swivel Glides

| Stretchers | Height | with Work Top | without Work Top |
|------------|--------|-----------------|------------------|
| No | 31" | T19W312454□00_0 | T19W312454-T00_0 |
| No | 36" | T19W362454□00_0 | T19W362454-T00_0 |
| Yes | 31" | T09W312454□00_0 | T09W312454-T00_0 |
| Yes | 36" | T09W362454□00_0 | T09W362454-T00_0 |

Wood Multi-Science Tables


Multi-Science Tables are available in both 31" high sitting and 36" high standing heights. They provide laboratory and lecture space for two or four students in teaching laboratories. Each table provides book/tote tray compartments in apron for each student, an electrical outlet in the apron rail and rod

sockets in table tops for ring support rods and/or upright rod assemblies (purchase rods separately). Work tops are 1" thick and overhang the apron on all four sides by 2". Legs are equipped with leg shoes and holddown angles for permanent placement.

Note: Work tops are available in three materials: Kemresin, Plastic Laminate and Acid Resistant Plastic Laminate. To specify top material, replace the box in part number with the top material designation letter.

Example: T13W312472R0010 specifies an Oak Two Student Table with a Kemresin work top.

General Science, Physics or Biology Student Tables



Two Student Table

24" Deep x 72" Long


Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate

Legs: Leg Shoes with Holddown Angles

Fixtures: (1) 0581-01 120 VAC Duplex

Fittings: (2) 0548-00 Rod Sockets

| Stretchers | Height | with Work Top Fixtures & Fittings | without Work Top Fixtures & Fittings |
|------------|--------|--------------------------------------|---|
| No | 31" | T13W312472□00_0 | T13W312472-T00_0 |
| No | 36" | T13W362472□00_0 | T13W362472-T00_0 |
| Yes | 31" | T03W312472□00_0 | T03W312472-T00_0 |
| Yes | 36" | T03W362472□00_0 | T03W362472-T00_0 |



Four Student Table

42" Deep x 72" Long


Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate

Legs: Leg Shoes with Holddown Angles

Fixtures: (2) 0581-01 120 VAC Duplex

Fittings: (4) 0548-00 Rod Sockets

| Stretchers | Height | with Work Top Fixtures & Fittings | without Work Top Fixtures & Fittings |
|------------|--------|--------------------------------------|---|
| No | 31" | T14W314272□00_0 | T14W314272-T00_0 |
| No | 36" | T14W364272□00_0 | T14W364272-T00_0 |
| Yes | 31" | T04W314272□00_0 | T04W314272-T00_0 |
| Yes | 36" | T04W364272□00_0 | T04W364272-T00_0 |



Two Student Table

24" Deep x 72" Long

Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate


Legs: Leg Shoes with Holddown Angles

Fixtures: (1) 0581-01 120 VAC Duplex

Fittings: (1) W-0263-00 Double Outlet Service Fitting

(2) 0548-00 Rod Sockets

| Stretchers | Height | with Work Top Fixtures & Fittings | without Work Top Fixtures & Fittings |
|------------|--------|--------------------------------------|---|
| No | 31" | T15W312472□00_0 | T15W312472-T00_0 |
| No | 36" | T15W362472□00_0 | T15W362472-T00_0 |
| Yes | 31" | T05W312472□00_0 | T05W312472-T00_0 |
| Yes | 36" | T05W362472□00_0 | T05W362472-T00_0 |



Four Student Table

42" Deep x 72" Long

Work Top: R – Black Kemresin
P – Black Plastic Laminate
D – Black Acid Resistant Plastic Laminate

Legs: Leg Shoes with Holddown Angles

Fixtures: (2) 0581-01 120 VAC Duplex

Fittings: (1) W-0272-00 Quadruple Outlet Service Fitting

(4) 0548-00 Rod Sockets

| Stretchers | Height | with Work Top Fixtures & Fittings | without Work Top Fixtures & Fittings |
|------------|--------|--------------------------------------|---|
| No | 31" | T16W314272□00_0 | T16W314272-T00_0 |
| No | 36" | T16W364272□00_0 | T16W364272-T00_0 |
| Yes | 31" | T06W314272□00_0 | T06W314272-T00_0 |
| Yes | 36" | T06W364272□00_0 | T06W364272-T00_0 |