

Wood Safety Cabinets

Safety Cabinet: a valuable addition for any laboratory to concentrate many safety equipment items in one location. Two half-depth fixed shelves and the lower cupboard provide storage space for fire blankets, fire extinguisher, first aid kit and supplies. The cabinet is equipped with an Emergency Shower Head Assembly and a Deck


Mounted Eyewash Assembly.

The **ADA Safety Cabinet** is the same except the Emergency Shower is fitted with a longer handle, and a Dual Purpose Eye Wash/Drench Hose is substituted for the Deck Mounted Eye Wash.

Both are furnished with base molding and corner clips.

Note: Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Safety Cabinet



G81W902236-_____ 90" High x 22" Deep x 36" Long


Cabinet: 90" H x 21¹/₄" D x 36" L with two fixed shelves and double cupboard

Eyewash: Deck Mounted with two Gentle Spray Outlet Heads, 1¹/₂" Dia. ABS Plastic Bowl and Stay-Open Ball Valve with Epoxy Coated Flag Handle

Shower: 10" Dia. ABS Plastic Shower Head with 1" I.P.S. Stay-Open Ball Valve with 29" long Pull Rod.

No service piping or electrical conduit included.

ADA Safety Cabinet



G84W902236-_____ 90" High x 22" Deep x 36" Long

Cabinet: 90" H x 21¹/₄" D x 36" L with two fixed shelves and double cupboard

Eyewash: **W-0934-00** Deck Mounted with two Gentle Spray Outlet Heads, Squeeze Valve with Locking Clip

Shower: 10" Dia. ABS Plastic Shower Head with 1" I.P.S. Stay-Open Ball Valve with 59¹/₂" long Pull Rod.

No service piping or electrical conduit included.


Sink Assemblies

These Sink Assemblies are designed for wash-up and general laboratory use. Each is furnished with epoxy resin worksurfaces, curbs, ledges and sinks. Drain grooves are supplied in the drainboards as illustrated. Filler panels, base molding and molding corner clips are

included to complete the sink assemblies.

Note: Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

38" Wall Sink Assembly



T68W363138-_____ 36" High x 31" Deep x 38" Long

Work Top: 1" Thick Black Kemresin
with 4" High Curb and Drain Grooves

Work Top Size: 31" x 38"

Cabinets: (1) G00W362236- Sink Cabinet

Sink: (1) **1003-00** Kemresin – 11" H x 15" W x 18" L

(1) **0482-BP** 1 1/2" I.P.S. Sink Outlet


Fittings: (1) **W-0340-0V** H & C Water w/VB Gooseneck

No service piping or electric conduit included.

Options:

T68W363138-T_____ Without Work Top, Sink & Fittings

50" Wall Sink Assembly



T68W363150-_____ 36" High x 31" Deep x 50" Long

Work Top: 1" Thick Black Kemresin
with 4" High Curb and Drain Grooves

Work Top Size: 31" x 50"

Cabinets: (1) G00W362248- Sink Cabinet

Sink: (1) **1005-00** Kemresin – 10" H x 15" W x 25" L

(1) **0482-BP** 1 1/2" I.P.S. Sink Outlet


Fittings: (1) **W-0340-0V** H & C Water w/VB Gooseneck

No service piping or electric conduit included.

Options:

T68W363150-T_____ Without Work Top, Sink & Fittings

Scrub Sink Assembly



T69W451925-_____ 45" High x 18 11/16" Deep x 25" Long

Ledge : 1 1/4" Thick Black Kemresin

Ledge Size: 3 1/4" x 25"

Height: Ledge 45" – Top of Sink 36"

Cabinet: (1) T69W231822-00_0 22 5/8" H x 17 11/16" D x 22" L

Sink: (1) **1037-00** Kemresin – 12 3/4" H x 14 1/4" W x 23 1/4" L

(1) **0482-BP** 1 1/2" I.P.S. Sink Outlet

Fittings: (1) **W-0355-00** Hot & Cold Water w/Swing Spout

No service piping or electric conduit included.