

STAINLESS STEEL BASIKBENCH



Engineered with the same strength, flexibility and cost-efficiency of the original BASIKBench, this new version of Kewaunee's classic can go where no BASIKBench has gone before.



Constructed of burnished, 304 stainless steel, this BASIKBench is perfect for the most demanding finish-free environments, including clean rooms, research and development laboratories and hospitals. And with plenty of customization options and accessories, it lives up to the Kewaunee name.

Engineered with our C-frame design, BASIKBench provides the option to design a custom workstation for maximum productivity in the workplace.

Available in 30" sit-down and 36" technical heights and 48", 60" and 72" widths, the stainless steel bench features options such as above-the-worksurface shelving, hanging storage and powered or non-powered shelf risers and worksurface backsplashes. BASIKBench's stand-alone design also allows for quick and easy reconfiguration of the workspace layout.



HYGIENE – Easy-to-clean stainless steel is ideal for hospitals, food processing, laboratories, clean rooms and any environment requiring finish-free products.



CORROSION-RESISTANCE – Rust-resistant 304 stainless steel is suitable for a wide range of atmospheric environments and many corrosive media.



DURABILITY – Stainless steel boasts superb heat resistance and an excellent strength-to-weight ratio, making it an intelligent, long-term investment. Moreover, a standard modesty panel and front apron provide even greater weight-bearing capability and support.



APPEARANCE – The burnished, easy-to-maintain surface of stainless steel graces the workspace with a modern, attractive appearance that won't fade with time.



www.kewaunee.com

2700 West Front Street
Statesville, NC 28677

Phone: 1-800-824-6626
Fax: 1-800-624-5536