

Floor Mounted Base Cabinets


Features:

Steel Base Cabinets


- Electrostatically applied, chemical resistant, powder coat finish available in colors to complement the laboratory environment.
- One piece drawer body-bottom, back, sides and inner drawer head, radiused all four sides.
- One-piece case bottom and horizontal front rail.
- Self-closing drawer slide with wrap around construction and nylon tired ball bearing rollers.
- Flush face corner construction.
- Interchangeable doors and drawers.
- Five-knuckle institutional type stainless steel hinge, 2³/₄" high.

Wood Base Cabinets

- Multiple coat, chemical resistant finish system formulated to highlight the natural grain patterns.
- Dovetailed drawer construction with nine-ply hardwood plywood drawer sides.
- 3/4" seven-ply hardwood plywood end panels, bottom and shelf, with solid hardwood facing.
- Epoxy-coated, self-closing drawer slides with nylon rollers and two stage positive hold-open stops.
- One-piece 1/4" hardboard back, dadoed into end panels and bottom and securely fastened.
- Institutional type five-knuckle offset door hinges rated for 150 pound loads.



Sitting Height Floor Mounted Base Cabinets



Specifications:


Steel Sitting Height Base Cabinets are 30" high and 22" deep. Wood cabinets are 30" high and a nominal 22" deep. Both available in lengths shown.

Open Cabinets have one full-depth adjustable shelf.


Cupboard cabinets have one full-depth adjustable shelf.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.


Open Base Cabinets




Steel Cabinet	Length	Wood Cabinet
B00M312215-0000	15"	B00W312215-00_0
B00M312218-0000	18"	B00W312218-00_0
B00M312221-0000	21"	B00W312221-00_0
B00M312224-0000	24"	B00W312224-00_0





Steel Cabinet	Length	Wood Cabinet
B00M312230-0000	30"	B00W312230-00_0
B00M312236-0000	36"	B00W312236-00_0
B00M312242-0000	42"	B00W312242-00_0
B00M312248-0000	48"	B00W312248-00_0





Cupboard Base Cabinets




Steel Cabinet	Length	Wood Cabinet
C03M312215L	15"	C03W312215-L
C03M312218L	18"	C03W312218-L
C03M312221L	21"	C03W312221-L
C03M312224L	24"	C03W312224-L




Steel Cabinet	Length	Wood Cabinet
C03M312215-	15"	C03W312215-
C03M312218-	18"	C03W312218-
C03M312221-	21"	C03W312221-
C03M312224-	24"	C03W312224-




Steel Cabinet	Length	Wood Cabinet
C03M312230-	30"	C03W312230-
C03M312236-	36"	C03W312236-
C03M312242-	42"	C03W312242-
C03M312248-	48"	C03W312248-




Drawer Base Cabinets



Steel Cabinet	Length	Wood Cabinet
D08M312215-	15"	D08W312215-_0
D08M312218-	18"	D08W312218-_0
D08M312221-	21"	D08W312221-_0
D08M312224-	24"	D08W312224-_0



Sitting Height Floor Mounted Base Cabinets




Specifications:



Steel Sitting Height Base Cabinets are 30" high and 22" deep. Wood cabinets are 30" high and a nominal 22" deep. Both available in lengths shown.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Drawer Base Cabinets





Steel Cabinet	Length	Wood Cabinet
D10M312215- ----	15"	D10W312215- 0 ----
D10M312218- ----	18"	D10W312218- 0 ----
D10M312221- ----	21"	D10W312221- 0 ----
D10M312224- ----	24"	D10W312224- 0 ----




Steel Cabinet	Length	Wood Cabinet
D30M312215- ----	15"	D30W312215- 0 ----
D30M312218- ----	18"	D30W312218- 0 ----
D30M312221- ----	21"	D30W312221- 0 ----
D30M312224- ----	24"	D30W312224- 0 ----



Steel Cabinet	Length	Wood Cabinet
D31M312230- ----	30"	D31W312230- 0 ----
D31M312236- ----	36"	D31W312236- 0 ----
D31M312242- ----	42"	D31W312242- 0 ----
D31M312248- ----	48"	D31W312248- 0 ----




Steel Cabinet	Length	Wood Cabinet
D81M312218- ----	18"	D81W312218- 0 ----
D81M312221- ----	21"	D81W312221- 0 ----
D81M312224- ----	24"	D81W312224- 0 ----




2 File Drawers with hangers

2 File Drawers with hangers




Steel Cabinet	Length
D83M312218- ----	18"
D83M312221- ----	21"
D83M312224- ----	24"

1 File Drawer with hangers




Steel Cabinet	Length	Wood Cabinet
D84M312218- ----	18"	D90W312218- 0 ----
D84M312221- ----	21"	D90W312221- 0 ----
D84M312224- ----	24"	D90W312224- 0 ----

Pullboard
1 File Drawer with hangers



Pullboard
1 File Drawer with hangers

Sitting Height Floor Mounted Base Cabinets




Specifications:


Steel Sitting Height Base Cabinets are 30" high and 22" deep. Wood cabinets are 30" high and a nominal 22" deep. Both available in lengths shown.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.


Drawer Base Cabinets





Steel Cabinet	Length	Wood Cabinet
D85M312215- -----	15"	D85W312215- 0 -----
D85M312218- -----	18"	D85W312218- 0 -----
D85M312221- -----	21"	D85W312221- 0 -----
D85M312224- -----	24"	D85W312224- 0 -----





Drawer & Cupboard Base Cabinets





Steel Cabinet	Length	Wood Cabinet
E40M312215L -----	15"	E40W312215L -----
E40M312218L -----	18"	E40W312218L -----
E40M312221L -----	21"	E40W312221L -----
E40M312224L -----	24"	E40W312224L -----




Steel Cabinet	Length	Wood Cabinet
E40M312215- -----	15"	E40W312215- -----
E40M312218- -----	18"	E40W312218- -----
E40M312221- -----	21"	E40W312221- -----
E40M312224- -----	24"	E40W312224- -----




Steel Cabinet	Length	Wood Cabinet
E44M312230L -----	30"	E44W312230L -----
E44M312236L -----	36"	E44W312236L -----
E44M312242L -----	42"	E44W312242L -----
E44M312248L -----	48"	E44W312248L -----



Steel Cabinet	Length	Wood Cabinet
E44M312230- -----	30"	E44W312230- -----
E44M312236- -----	36"	E44W312236- -----
E44M312242- -----	42"	E44W312242- -----
E44M312248- -----	48"	E44W312248- -----



ADA Height Floor Mounted Base Cabinets



Specifications:


Steel ADA Height Base Cabinets are 33" high and 22" deep. Wood cabinets are 32½" high and a nominal 22" deep. Both available in lengths shown.

Open Cabinets have one full-depth adjustable shelf.



Cupboard cabinets have one full-depth adjustable shelf.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Cupboard Base Cabinets




Steel Cabinet	Length	Wood Cabinet
C03M342215L -----	15"	C03W342215L -----
C03M342218L -----	18"	C03W342218L -----
C03M342221L -----	21"	C03W342221L -----
C03M342224L -----	24"	C03W342224L -----




Steel Cabinet	Length	Wood Cabinet
C03M342215- -----	15"	C03W342215- -----
C03M342218- -----	18"	C03W342218- -----
C03M342221- -----	21"	C03W342221- -----
C03M342224- -----	24"	C03W342224- -----




Steel Cabinet	Length	Wood Cabinet
C03M342230- -----	30"	C03W342230- -----
C03M342236- -----	36"	C03W342236- -----
C03M342242- -----	42"	C03W342242- -----
C03M342248- -----	48"	C03W342248- -----




Drawer Base Cabinets





Steel Cabinet	Length	Wood Cabinet
D30M342215- -----	15"	D30W342215- 0 -----
D30M342218- -----	18"	D30W342218- 0 -----
D30M342221- -----	21"	D30W342221- 0 -----
D30M342224- -----	24"	D30W342224- 0 -----




Drawer & Cupboard Base Cabinets




Steel Cabinet	Length	Wood Cabinet
E40M342215L -----	15"	E40W342215L -----
E40M342218L -----	18"	E40W342218L -----
E40M342221L -----	21"	E40W342221L -----
E40M342224L -----	24"	E40W342224L -----



Steel Cabinet	Length	Wood Cabinet
E40M342215- -----	15"	E40W342215- -----
E40M342218- -----	18"	E40W342218- -----
E40M342221- -----	21"	E40W342221- -----
E40M342224- -----	24"	E40W342224- -----



ADA Height Floor Mounted Base Cabinets



Specifications:


Steel ADA Height Base Cabinets are 33" high and 22" deep. Wood cabinets are 32½" high and a nominal 22" deep. Both available in lengths shown.

Open Cabinets have one full-depth adjustable shelf.



Cupboard cabinets have one full-depth adjustable shelf.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Drawer & Cupboard Base Cabinets





Steel Cabinet	Length	Wood Cabinet
E40M342230-	30"	E40W342230-
E40M342236-	36"	E40W342236-
E40M342242-	42"	E40W342242-
E40M342248-	48"	E40W342248-





Steel Cabinet	Length	Wood Cabinet
E41M342230-	30"	E41W342230-
E41M342236-	36"	E41W342236-
E41M342242-	42"	E41W342242-
E41M342248-	48"	E41W342248-





Steel Cabinet	Length	Wood Cabinet
E43M342230L	30"	E43W342230L
E43M342236L	36"	E43W342236L
E43M342242L	42"	E43W342242L
E43M342248L	48"	E43W342248L




Steel Cabinet	Length	Wood Cabinet
E43M342230-	30"	E43W342230-
E43M342236-	36"	E43W342236-
E43M342242-	42"	E43W342242-
E43M342248-	48"	E43W342248-




Steel Cabinet	Length	Wood Cabinet
E45M342230L	30"	E45W342230L
E45M342236L	36"	E45W342236L
E45M342242L	42"	E45W342242L
E45M342248L	48"	E45W342248L



Steel Cabinet	Length	Wood Cabinet
E45M342230-	30"	E45W342230-
E45M342236-	36"	E45W342236-
E45M342242-	42"	E45W342242-
E45M342248-	48"	E45W342248-



ADA Height Floor Mounted Base Cabinets




Specifications:



Steel ADA Height Base Cabinets are 33" high and 22" deep. Wood cabinets are 32½" high and a nominal 22" deep. Both available in lengths shown.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Sink Base Cabinets





Steel Cabinet	Length	Wood Cabinet
G00M342215L -----	15"	G00W342215L -----
G00M342218L -----	18"	G00W342218L -----
G00M342221L -----	21"	G00W342221L -----
G00M342224L -----	24"	G00W342224L -----




Steel Cabinet	Length	Wood Cabinet
G00M342215- -----	15"	G00W342215- -----
G00M342218- -----	18"	G00W342218- -----
G00M342221- -----	21"	G00W342221- -----
G00M342224- -----	24"	G00W342224- -----




Steel Cabinet	Length	Wood Cabinet
G00M342230- -----	30"	G00W342230- -----
G00M342236- -----	36"	G00W342236- -----
G00M342242- -----	42"	G00W342242- -----
G00M342248- -----	48"	G00W342248- -----



Steel Cabinet	Length	Wood Cabinet
G03M342230- -----	30"	G03W342230- -----
G03M342236- -----	36"	G03W342236- -----
G03M342242- -----	42"	G03W342242- -----
G03M342248- -----	48"	G03W342248- -----



Standing Height Floor Mounted Cabinets



Specifications:


Steel Standing Height Base Cabinets are 35 3/4" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.

Open Cabinets have two full-depth adjustable shelves.



Cupboard cabinets have one full-depth adjustable shelf.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.


Open Base Cabinets




Steel Cabinet	Length	Wood Cabinet
B00M372215-0000	15"	B00W362215-00_0
B00M372218-0000	18"	B00W362218-00_0
B00M372221-0000	21"	B00W362221-00_0
B00M372224-0000	24"	B00W362224-00_0





Steel Cabinet	Length	Wood Cabinet
B00M372230-0000	30"	B00W362230-00_0
B00M372236-0000	36"	B00W362236-00_0
B00M372242-0000	42"	B00W362242-00_0
B00M372248-0000	48"	B00W362248-00_0




Cupboard Base Cabinets




Steel Cabinet	Length	Wood Cabinet
C03M372215L	15"	C03W362215L
C03M372218L	18"	C03W362218L
C03M372221L	21"	C03W362221L
C03M372224L	24"	C03W362224L




Steel Cabinet	Length	Wood Cabinet
C03M372215-	15"	C03W362215-
C03M372218-	18"	C03W362218-
C03M372221-	21"	C03W362221-
C03M372224-	24"	C03W362224-




Steel Cabinet	Length	Wood Cabinet
C03M372230-	30"	C03W362230-
C03M372236-	36"	C03W362236-
C03M372242-	42"	C03W362242-
C03M372248-	48"	C03W362248-




Drawer Base Cabinets



Steel Cabinet	Length	Wood Cabinet
D30M372215-	15"	D30W362215-_0
D30M372218-	18"	D30W362218-_0
D30M372221-	21"	D30W362221-_0
D30M372224-	24"	D30W362224-_0



Standing Height Floor Mounted Cabinets




Specifications:



Steel Standing Height Base Cabinets are 35³/₄" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Drawer Base Cabinets





Steel Cabinet	Length	Wood Cabinet
D30M372230- ----	30"	D30W362230- 0
D30M372236- ----	36"	D30W362236- 0
D30M372242- ----	42"	D30W362242- 0
D30M372248- ----	48"	D30W362248- 0





Steel Cabinet	Length	Wood Cabinet
D31M372230- ----	30"	D31W362230- 0
D31M372236- ----	36"	D31W362236- 0
D31M372242- ----	42"	D31W362242- 0
D31M372248- ----	48"	D31W362248- 0





Steel Cabinet	Length	Wood Cabinet
D38M372215- ----	15"	D38W362215- 0
D38M372218- ----	18"	D38W362218- 0
D38M372221- ----	21"	D38W362221- 0
D38M372224- ----	24"	D38W362224- 0




Steel Cabinet	Length	Wood Cabinet
D39M372230- ----	30"	D39W362230- 0
D39M372236- ----	36"	D39W362236- 0
D39M372242- ----	42"	D39W362242- 0
D39M372248- ----	48"	D39W362248- 0




Steel Cabinet	Length	Wood Cabinet
D50M372215- ----	15"	D50W362215- 0
D50M372218- ----	18"	D50W362218- 0
D50M372221- ----	21"	D50W362221- 0
D50M372224- ----	24"	D50W362224- 0



Steel Cabinet	Length	Wood Cabinet
D51M372230- ----	30"	D51W362230- 0
D51M372236- ----	36"	D51W362236- 0
D51M372242- ----	42"	D51W362242- 0
D51M372248- ----	48"	D51W362248- 0



Standing Height Floor Mounted Cabinets




Specifications:

Steel Standing Height Base Cabinets are 35 3/4" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.


See page 95 for designating cabinet style, door and drawer head material, and hardware selection.

Drawer Base Cabinets




2 File Drawers with hangers

Steel Cabinet	Length	Wood Cabinet
D81M372215- -----	15"	D81W362215- 0 -----
D81M372218- -----	18"	D81W362218- 0 -----
D81M372221- -----	21"	D81W362221- 0 -----
D81M372224- -----	24"	D81W362224- 0 -----




2 File Drawers with hangers




Pullboard
1 File Drawer with hangers

Steel Cabinet	Length	Wood Cabinet
D84M372215- -----	15"	D84W362215- 0 -----
D84M372218- -----	18"	D84W362218- 0 -----
D84M372221- -----	21"	D84W362221- 0 -----
D84M372224- -----	24"	D84W362224- 0 -----





Pullboard
1 File Drawer with hangers


Drawer & Cupboard Base Cabinets




Steel Cabinet	Length	Wood Cabinet
E40M372215L -----	15"	E40W362215L -----
E40M372218L -----	18"	E40W362218L -----
E40M372221L -----	21"	E40W362221L -----
E40M372224L -----	24"	E40W362224L -----




Steel Cabinet	Length	Wood Cabinet
E40M372215- -----	15"	E40W362215- -----
E40M372218- -----	18"	E40W362218- -----
E40M372221- -----	21"	E40W362221- -----
E40M372224- -----	24"	E40W362224- -----



Steel Cabinet	Length	Wood Cabinet
E40M372230- -----	30"	E40W362230- -----
E40M372236- -----	36"	E40W362236- -----
E40M372242- -----	42"	E40W362242- -----
E40M372248- -----	48"	E40W362248- -----



Standing Height Floor Mounted Cabinets




Specifications:

Steel Standing Height Base Cabinets are 35³/₄" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.



Cupboards have one full-depth adjustable shelf.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Drawer & Cupboard Base Cabinets





Steel Cabinet	Length	Wood Cabinet
E41M372230-	30"	E41W362230-
E41M372236-	36"	E41W362236-
E41M372242-	42"	E41W362242-
E41M372248-	48"	E41W362248-





Steel Cabinet	Length	Wood Cabinet
E43M372230L	30"	E43W362230L
E43M372236L	36"	E43W362236L
E43M372242L	42"	E43W362242L
E43M372248L	48"	E43W362248L




Steel Cabinet	Length	Wood Cabinet
E43M372230-	30"	E43W362230-
E43M372236-	36"	E43W362236-
E43M372242-	42"	E43W362242-
E43M372248-	48"	E43W362248-




Steel Cabinet	Length	Wood Cabinet
E45M372230L	30"	E45W362230L
E45M372236L	36"	E45W362236L
E45M372242L	42"	E45W362242L
E45M372248L	48"	E45W362248L



Steel Cabinet	Length	Wood Cabinet
E45M372230-	30"	E45W362230-
E45M372236-	36"	E45W362236-
E45M372242-	42"	E45W362242-
E45M372248-	48"	E45W362248-



Standing Height Floor Mounted Cabinets




Specifications:

Steel Standing Height Base Cabinets are 35 3/4" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.



Cupboards have one full-depth adjustable shelf.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Drawer & Cupboard Base Cabinets





Steel Cabinet	Length	Wood Cabinet
E50M372215L	15"	E52W362215L
E50M372218L	18"	E52W362218L
E50M372221L	21"	E52W362221L
E50M372224L	24"	E52W362224L





Steel Cabinet	Length	Wood Cabinet
E50M372215-	15"	E52W362215-
E50M372218-	18"	E52W362218-
E50M372221-	21"	E52W362221-
E50M372224-	24"	E52W362224-




Steel Cabinet	Length	Wood Cabinet
E51M372230-	30"	E53W362230-
E51M372236-	36"	E53W362236-
E51M372242-	42"	E53W362242-
E51M372248-	48"	E53W362248-




Steel Cabinet	Length	Wood Cabinet
E63M372230L	30"	E63W362230L
E63M372236L	36"	E63W362236L
E63M372242L	42"	E63W362242L
E63M372248L	48"	E63W362248L



Steel Cabinet	Length	Wood Cabinet
E63M372230-	30"	E63W362230-
E63M372236-	36"	E63W362236-
E63M372242-	42"	E63W362242-
E63M372248-	48"	E63W362248-



Standing Height Floor Mounted Cabinets




Specifications:

Steel Standing Height Base Cabinets are 35³/₄" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.



Cupboards have one full-depth adjustable shelf.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.


Drawer & Cupboard Base Cabinets




Steel Cabinet	Length	Wood Cabinet
E65M372230L	30"	E65W362230L
E65M372236L	36"	E65W362236L
E65M372242L	42"	E65W362242L
E65M372248L	48"	E65W362248L





Steel Cabinet	Length	Wood Cabinet
E65M372230-	30"	E65W362230-
E65M372236-	36"	E65W362236-
E65M372242-	42"	E65W362242-
E65M372248-	48"	E65W362248-





Sink Base Cabinets




Steel Cabinet	Length	Wood Cabinet
G00M372215L	15"	G00W362215L
G00M372218L	18"	G00W362218L
G00M372221L	21"	G00W362221L
G00M372224L	24"	G00W362224L




Steel Cabinet	Length	Wood Cabinet
G00M372215-	15"	G00W362215-
G00M372218-	18"	G00W362218-
G00M372221-	21"	G00W362221-
G00M372224-	24"	G00W362224-



Steel Cabinet	Length	Wood Cabinet
G00M372230-	30"	G00W362230-
G00M372236-	36"	G00W362236-
G00M372242-	42"	G00W362242-
G00M372248-	48"	G00W362248-



Standing Height Floor Mounted Cabinets




Specifications:

Steel Standing Height Base Cabinets are 35 3/4" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.



Sink cabinets are provided with sink supports.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.


Sink Base Cabinets




Steel Cabinet	Length	Wood Cabinet
G00M372260-	60"	G01W362260-





Steel Cabinet	Length	Wood Cabinet
G04M372252-	52"	G04W362252-
G04M372258-	58"	G04W362258-
G04M372260-	60"	G04W362260-





Fume Hood Base Cabinets





Steel Cabinet	Length	Wood Cabinet
G08M372215L	15"	G08W362215L
G08M372218L	18"	G08W362218L
G08M372221L	21"	G08W362221L
G08M372224L	24"	G08W362224L




Steel Cabinet	Length	Wood Cabinet
G08M372215-	15"	G08W362215-
G08M372218-	18"	G08W362218-
G08M372221-	21"	G08W362221-
G08M372224-	24"	G08W362224-




Steel Cabinet	Length	Wood Cabinet
G08M372230-	30"	G08W362230-
G08M372236-	36"	G08W362236-
G08M372242-	42"	G08W362242-
G08M372248-	48"	G08W362248-



Steel Cabinet	Length	Wood Cabinet
G08M372260-	60"	G08W362260-



Standing Height Floor Mounted Cabinets



Specifications:


Steel Standing Height Base Cabinets are 35³/₄" high and 22" deep. Wood cabinets are 35" high and a nominal 22" deep. Both available in lengths shown.

Fume Hood Acid Storage Base Cabinets are manufactured with a


one-piece polyethylene corrosion resistant lining and are provided with a PVC vent.

See page 95 for designating cabinet style, door and drawer head material, and hardware selection.



Fume Hood Acid Storage Base Cabinets





Steel Cabinet	Length	Wood Cabinet
G80M372224L _____	24"	G80W362224L _____




Steel Cabinet	Length	Wood Cabinet
G80M372224- _____	24"	G80W362224- _____




Steel Cabinet	Length	Wood Cabinet
G80M372230- _____	30"	G80W362230- _____
G80M372236- _____	36"	G80W362236- _____
G80M372248- _____	48"	G80W362248- _____



Steel Cabinet	Length	Wood Cabinet
G80M372260- _____	60"	G80W362260- _____




Diagonal Corner Base Cabinets




Sink Cabinet


Steel Cabinet	Length	Wood Cabinet
G14M372214- _____	12" - 32"	G14W362214- _____




Sink Cabinet



Steel Cabinet	Length	Wood Cabinet
G16M372214- _____	12" - 32"	G16W362214- _____



Standing Height Floor Mounted Cabinets




Specifications:



Steel Solvent Storage Base Cabinets, are 35 3/4" high, and 22" deep. Wood Solvent Storage Base Cabinets, are 35" high, and a nominal 22" deep. Lengths as shown. They are specifically designed for the storage of flammable and combustible liquids.

These UL listed Solvent Storage Cabinets are provided with self-closing, self-latching doors, and meet UFC, OSHA and NFPA No.30-1993 requirements. Steel available in Style 1 only, Wood available in Styles 1 & 5 only.



Solvent Storage Base Cabinets





Steel Cabinet	Length	Wood Cabinet
G63M3724L0000	24"	G63W362224L_0_0





Steel Cabinet	Length	Wood Cabinet
G63M3724-0000	24"	G63W362224-_0_0





Steel Cabinet	Length	Wood Cabinet
G63M3730-0000	30"	G63W362230-_0_0
G63M3736-0000	36"	G63W362236-_0_0
G63M3748-0000	48"	G63W362248-_0_0



Steel Cabinet	Length	Wood Cabinet
G68M3724L0000	24"	G68W362224L_0_0



Steel Cabinet	Length	Wood Cabinet
G68M3724-0000	24"	G68W362224-_0_0



Steel Cabinet	Length	Wood Cabinet
G68M3730-0000	30"	G68W362230-_0_0
G68M3736-0000	36"	G68W362236-_0_0
G68M3748-0000	48"	G68W362248-_0_0

