


Upper & Full Height Storage Cabinets


Features:

- Both Steel and Wood Upper Cabinets are available.
- Upper Cabinets are available 24" and 30" high, and 12³/₄" and 16" deep.
- Same Upper Cabinet may be mounted on wall or on Alpha System.
- Upper Cabinets may span two Structural Modules.
- Multiple Upper Cabinets may be mounted on a single Structural Module.
- Both Steel and Wood Full Height Storage Cabinets are available.
- Full Height Storage Cabinets match height of Upper Carrier Structural Modules.
- Full Height Storage Cabinets available 12³/₄", 16" or 22" deep.



Upper & Full Height Storage Cabinets

Upper Cabinet Hanging System




With the use of an Upper Hanging Bracket and a pair of Hanging Rails, Kewaunee's standard Wall Cabinets can be mounted on Alpha System's Structural Module. Two Hanging Rails are snapped into the Upper Carrier upright slots. The Upper Hanging Bracket bolts to the top of the upper cabinet which then hangs on the upper Hanging Rail. The lower Hanging Rail allows the bottom of the cabinet to be securely fastened to the module. This system allows the Upper Cabinet to be moved both vertically and horizontally on the Structural Module.

Hanging Rails and Upper Hanging Brackets, are finished to match the Structural Module and furnished with all required hardware.

Hanging Rails are furnished in pairs.


Order the appropriate Upper Hanging Bracket in the same length as the Upper Cabinet. The Hanging Rails must be the same length as the Structural Module.

Upper Hanging Bracket




Length	Part Number
18"	F48M0418-00
24"	F48M0424-00
30"	F48M0430-00
36"	F48M0436-00
42"	F48M0442-00
48"	F48M0448-00

Hanging Rails



Length	Part Number	
24"	F49M0224-00	1 pair
30"	F49M0230-00	1 pair
36"	F49M0236-00	1 pair
42"	F49M0242-00	1 pair
48"	F49M0248-00	1 pair
54"	F49M0254-00	1 pair
60"	F49M0260-00	1 pair
66"	F49M0266-00	1 pair
72"	F49M0272-00	1 pair

Upper Cabinets – 24" High



Specifications:


These Upper Cabinets are 24" high and available in two depths. Lengths as shown. They are provided with one adjustable shelf, and are available in either wood or steel.

Sliding Plate Glass doors are mounted on friction-free roller glides located at the bottom of the door.



See page 95 for designating cabinet material on wood wall cabinets.

Hanging system must be ordered separately.


Open Storage Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W00M241315-0000	15"	W00W241215-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M241318-0000	18"	W00W241218-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M241321-0000	21"	W00W241221-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M241324-0000	24"	W00W241224-00_0	12 ³ / ₄ "
16"	W00M241615-0000	15"	W00W241615-00_0	16 ³ / ₄ "
16"	W00M241618-0000	18"	W00W241618-00_0	16 ³ / ₄ "
16"	W00M241621-0000	21"	W00W241621-00_0	16 ³ / ₄ "
16"	W00M241624-0000	24"	W00W241624-00_0	16 ³ / ₄ "




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W00M241330-0000	30"	W00W241230-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M241336-0000	36"	W00W241236-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M241342-0000	42"	W00W241242-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M241348-0000	48"	W00W241248-00_0	12 ³ / ₄ "
16"	W00M241630-0000	30"	W00W241630-00_0	16 ³ / ₄ "
16"	W00M241636-0000	36"	W00W241636-00_0	16 ³ / ₄ "
16"	W00M241642-0000	42"	W00W241642-00_0	16 ³ / ₄ "
16"	W00M241648-0000	48"	W00W241648-00_0	16 ³ / ₄ "




Sliding Plate Glass Door Cabinets



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W10M241330-0000	30"	W10W241230-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W10M241336-0000	36"	W10W241236-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W10M241342-0000	42"	W10W241242-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W10M241348-0000	48"	W10W241248-00_0	12 ³ / ₄ "
16"	W10M241630-0000	30"	W10W241630-00_0	16 ³ / ₄ "
16"	W10M241636-0000	36"	W10W241636-00_0	16 ³ / ₄ "
16"	W10M241642-0000	42"	W10W241642-00_0	16 ³ / ₄ "
16"	W10M241648-0000	48"	W10W241648-00_0	16 ³ / ₄ "




Sliding Plate Glass Door Cabinets with Smoked Glass Doors



Depth	Steel Cabinet	Length
12 ³ / ₄ "	W30M241330-0000	30"
12 ³ / ₄ "	W30M241336-0000	36"
12 ³ / ₄ "	W30M241342-0000	42"
12 ³ / ₄ "	W30M241348-0000	48"
16"	W30M241630-0000	30"
16"	W30M241636-0000	36"
16"	W30M241642-0000	42"
16"	W30M241648-0000	48"

Upper Cabinets – 24" High



Specifications:


These Upper Cabinets are 24" high and available in two depths. Lengths as shown. They are provided with one adjustable shelf, and are available in either wood or steel.

Sliding Glazed and Panel Doors are suspended by roller from a track at the top of the cabinet.


See page 95 for designating cabinet door style, door and cabinet material, and pull and hinge selection.

Hanging system must be ordered separately.


Sliding Panel Door Cabinets





Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W12M241330-__0	15"	W12W241230-__0_0	12 ³ / ₄ "
12 ³ / ₄ "	W12M241336-__0	18"	W12W241236-__0_0	12 ³ / ₄ "
12 ³ / ₄ "	W12M241342-__0	21"	W12W241242-__0_0	12 ³ / ₄ "
12 ³ / ₄ "	W12M241348-__0	24"	W12W241248-__0_0	12 ³ / ₄ "
16"	W12M241630-__0	15"	W12W241630-__0_0	16 ³ / ₄ "
16"	W12M241636-__0	18"	W12W241636-__0_0	16 ³ / ₄ "
16"	W12M241642-__0	21"	W12W241642-__0_0	16 ³ / ₄ "
16"	W12M241648-__0	24"	W12W241648-__0_0	16 ³ / ₄ "





Swinging Glazed Door Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W20M241315L	15"	W20W241215L	12 ³ / ₄ "
12 ³ / ₄ "	W20M241318L	18"	W20W241218L	12 ³ / ₄ "
12 ³ / ₄ "	W20M241321L	21"	W20W241221L	12 ³ / ₄ "
12 ³ / ₄ "	W20M241324L	24"	W20W241224L	12 ³ / ₄ "
16"	W20M241615L	15"	W20W241615L	16 ³ / ₄ "
16"	W20M241618L	18"	W20W241618L	16 ³ / ₄ "
16"	W20M241621L	21"	W20W241621L	16 ³ / ₄ "
16"	W20M241624L	24"	W20W241624L	16 ³ / ₄ "




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W20M241315-	15"	W20W241215-	12 ³ / ₄ "
12 ³ / ₄ "	W20M241318-	18"	W20W241218-	12 ³ / ₄ "
12 ³ / ₄ "	W20M241321-	21"	W20W241221-	12 ³ / ₄ "
12 ³ / ₄ "	W20M241324-	24"	W20W241224-	12 ³ / ₄ "
16"	W20M241615-	15"	W20W241615-	16 ³ / ₄ "
16"	W20M241618-	18"	W20W241618-	16 ³ / ₄ "
16"	W20M241621-	21"	W20W241621-	16 ³ / ₄ "
16"	W20M241624-	24"	W20W241624-	16 ³ / ₄ "



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W20M241330-	30"	W20W241230-	12 ³ / ₄ "
12 ³ / ₄ "	W20M241336-	36"	W20W241236-	12 ³ / ₄ "
12 ³ / ₄ "	W20M241342-	42"	W20W241242-	12 ³ / ₄ "
12 ³ / ₄ "	W20M241348-	48"	W20W241248-	12 ³ / ₄ "
16"	W20M241630-	30"	W20W241630-	16 ³ / ₄ "
16"	W20M241636-	36"	W20W241636-	16 ³ / ₄ "
16"	W20M241642-	42"	W20W241642-	16 ³ / ₄ "
16"	W20M241648-	48"	W20W241648-	16 ³ / ₄ "



Upper Cabinets – 24" High




Specifications:

These Upper Cabinets are 24" high and available in two depths. Lengths as shown. They are provided with one adjustable shelf, and are available in either wood or steel.



See page 95 for designating cabinet door style, door and cabinet material, and pull and hinge selection.

Hanging system must be ordered separately.



Swinging Panel Door Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W25M241315L	15"	W25W241215L	12 ³ / ₄ "
12 ³ / ₄ "	W25M241318L	18"	W25W241218L	12 ³ / ₄ "
12 ³ / ₄ "	W25M241321L	21"	W25W241221L	12 ³ / ₄ "
12 ³ / ₄ "	W25M241324L	24"	W25W241224L	12 ³ / ₄ "
16"	W25M241615L	15"	W25W241615L	16 ³ / ₄ "
16"	W25M241618L	18"	W25W241618L	16 ³ / ₄ "
16"	W25M241621L	21"	W25W241621L	16 ³ / ₄ "
16"	W25M241624L	24"	W25W241624L	16 ³ / ₄ "




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W25M241315-	15"	W25W241215-	12 ³ / ₄ "
12 ³ / ₄ "	W25M241318-	18"	W25W241218-	12 ³ / ₄ "
12 ³ / ₄ "	W25M241321-	21"	W25W241221-	12 ³ / ₄ "
12 ³ / ₄ "	W25M241324-	24"	W25W241224-	12 ³ / ₄ "
16"	W25M241615-	15"	W25W241615-	16 ³ / ₄ "
16"	W25M241618-	18"	W25W241618-	16 ³ / ₄ "
16"	W25M241621-	21"	W25W241621-	16 ³ / ₄ "
16"	W25M241624-	24"	W25W241624-	16 ³ / ₄ "



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W25M241330-	30"	W25W241230-	12 ³ / ₄ "
12 ³ / ₄ "	W25M241336-	36"	W25W241236-	12 ³ / ₄ "
12 ³ / ₄ "	W25M241342-	42"	W25W241242-	12 ³ / ₄ "
12 ³ / ₄ "	W25M241348-	48"	W25W241248-	12 ³ / ₄ "
16"	W25M241630-	30"	W25W241630-	16 ³ / ₄ "
16"	W25M241636-	36"	W25W241636-	16 ³ / ₄ "
16"	W25M241642-	42"	W25W241642-	16 ³ / ₄ "
16"	W25M241648-	48"	W25W241648-	16 ³ / ₄ "



Upper Cabinets – 30" High



Specifications:


These Upper Cabinets are 30" high and available in two depths. Lengths as shown. They are provided with two adjustable shelves, and are available in either wood or steel.

Sliding Plate Glass doors are mounted on friction-free roller glides located at the bottom of the door.



See page 95 for designating cabinet material on wood wall cabinets.

Hanging system must be ordered separately.

Open Storage Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W00M301315-0000	15"	W00W301215-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M301318-0000	18"	W00W301218-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M301321-0000	21"	W00W301221-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M301324-0000	24"	W00W301224-00_0	12 ³ / ₄ "
16"	W00M301615-0000	15"	W00W301615-00_0	16 ³ / ₄ "
16"	W00M301618-0000	18"	W00W301618-00_0	16 ³ / ₄ "
16"	W00M301621-0000	21"	W00W301621-00_0	16 ³ / ₄ "
16"	W00M301624-0000	24"	W00W301624-00_0	16 ³ / ₄ "




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W00M301330-0000	30"	W00W301230-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M301336-0000	36"	W00W301236-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M301342-0000	42"	W00W301242-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W00M301348-0000	48"	W00W301248-00_0	12 ³ / ₄ "
16"	W00M301630-0000	30"	W00W301630-00_0	16 ³ / ₄ "
16"	W00M301636-0000	36"	W00W301636-00_0	16 ³ / ₄ "
16"	W00M301642-0000	42"	W00W301642-00_0	16 ³ / ₄ "
16"	W00M301648-0000	48"	W00W301648-00_0	16 ³ / ₄ "




Sliding Plate Glass Door Cabinets



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W10M301330-0000	30"	W10W301230-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W10M301336-0000	36"	W10W301236-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W10M301342-0000	42"	W10W301242-00_0	12 ³ / ₄ "
12 ³ / ₄ "	W10M301348-0000	48"	W10W301248-00_0	12 ³ / ₄ "
16"	W10M301630-0000	30"	W10W301630-00_0	16 ³ / ₄ "
16"	W10M301636-0000	36"	W10W301636-00_0	16 ³ / ₄ "
16"	W10M301642-0000	42"	W10W301642-00_0	16 ³ / ₄ "
16"	W10M301648-0000	48"	W10W301648-00_0	16 ³ / ₄ "




Sliding Plate Glass Door Cabinets with Smoked Glass Doors



Depth	Steel Cabinet	Length
12 ³ / ₄ "	W30M301330-0000	30"
12 ³ / ₄ "	W30M301336-0000	36"
12 ³ / ₄ "	W30M301342-0000	42"
12 ³ / ₄ "	W30M301348-0000	48"
16"	W30M301630-0000	30"
16"	W30M301636-0000	36"
16"	W30M301642-0000	42"
16"	W30M301648-0000	48"

Upper Cabinets – 30" High



Specifications:


These Upper Cabinets are 30" high and available in two depths. Lengths as shown. They are provided with two adjustable shelves, and are available in either wood or steel.

Sliding Glazed and Panel Doors are suspended by roller from a track at the top of the cabinet.


See page 95 for designating cabinet door style, door and cabinet material, and pull and hinge selection.

Hanging system must be ordered separately.


Sliding Panel Door Cabinets





Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W12M301330-___0	15"	W12W301230-__0_0	12 ³ / ₄ "
12 ³ / ₄ "	W12M301336-___0	18"	W12W301236-__0_0	12 ³ / ₄ "
12 ³ / ₄ "	W12M301342-___0	21"	W12W301242-__0_0	12 ³ / ₄ "
12 ³ / ₄ "	W12M301348-___0	24"	W12W301248-__0_0	12 ³ / ₄ "
16"	W12M301630-___0	15"	W12W301630-__0_0	16 ³ / ₄ "
16"	W12M301636-___0	18"	W12W301636-__0_0	16 ³ / ₄ "
16"	W12M301642-___0	21"	W12W301642-__0_0	16 ³ / ₄ "
16"	W12M301648-___0	24"	W12W301648-__0_0	16 ³ / ₄ "





Swinging Glazed Door Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W20M301315L	15"	W20W301215L	12 ³ / ₄ "
12 ³ / ₄ "	W20M301318L	18"	W20W301218L	12 ³ / ₄ "
12 ³ / ₄ "	W20M301321L	21"	W20W301221L	12 ³ / ₄ "
12 ³ / ₄ "	W20M301324L	24"	W20W301224L	12 ³ / ₄ "
16"	W20M301615L	15"	W20W301615L	16 ³ / ₄ "
16"	W20M301618L	18"	W20W301618L	16 ³ / ₄ "
16"	W20M301621L	21"	W20W301621L	16 ³ / ₄ "
16"	W20M301624L	24"	W20W301624L	16 ³ / ₄ "




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W20M301315-	15"	W20W301215-	12 ³ / ₄ "
12 ³ / ₄ "	W20M301318-	18"	W20W301218-	12 ³ / ₄ "
12 ³ / ₄ "	W20M301321-	21"	W20W301221-	12 ³ / ₄ "
12 ³ / ₄ "	W20M301324-	24"	W20W301224-	12 ³ / ₄ "
16"	W20M301615-	15"	W20W301615-	16 ³ / ₄ "
16"	W20M301618-	18"	W20W301618-	16 ³ / ₄ "
16"	W20M301621-	21"	W20W301621-	16 ³ / ₄ "
16"	W20M301624-	24"	W20W301624-	16 ³ / ₄ "



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W20M301330-___	30"	W20W301230-___	12 ³ / ₄ "
12 ³ / ₄ "	W20M301336-___	36"	W20W301236-___	12 ³ / ₄ "
12 ³ / ₄ "	W20M301342-___	42"	W20W301242-___	12 ³ / ₄ "
12 ³ / ₄ "	W20M301348-___	48"	W20W301248-___	12 ³ / ₄ "
16"	W20M301630-___	30"	W20W301630-___	16 ³ / ₄ "
16"	W20M301636-___	36"	W20W301636-___	16 ³ / ₄ "
16"	W20M301642-___	42"	W20W301642-___	16 ³ / ₄ "
16"	W20M301648-___	48"	W20W301648-___	16 ³ / ₄ "



Upper Cabinets – 30" High




Specifications:

These Upper Cabinets are 30" high and available in two depths. Lengths as shown. They are provided with two adjustable shelves, and are available in either wood or steel.



See page 95 for designating cabinet door style, door and cabinet material, and pull and hinge selection.

Hanging system must be ordered separately.



Swinging Panel Door Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W25M301315L	15"	W25W301215L	12 ³ / ₄ "
12 ³ / ₄ "	W25M301318L	18"	W25W301218L	12 ³ / ₄ "
12 ³ / ₄ "	W25M301321L	21"	W25W301221L	12 ³ / ₄ "
12 ³ / ₄ "	W25M301324L	24"	W25W301224L	12 ³ / ₄ "
16"	W25M301615L	15"	W25W301615L	16 ³ / ₄ "
16"	W25M301618L	18"	W25W301618L	16 ³ / ₄ "
16"	W25M301621L	21"	W25W301621L	16 ³ / ₄ "
16"	W25M301624L	24"	W25W301624L	16 ³ / ₄ "




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W25M301315-	15"	W25W301215-	12 ³ / ₄ "
12 ³ / ₄ "	W25M301318-	18"	W25W301218-	12 ³ / ₄ "
12 ³ / ₄ "	W25M301321-	21"	W25W301221-	12 ³ / ₄ "
12 ³ / ₄ "	W25M301324-	24"	W25W301224-	12 ³ / ₄ "
16"	W25M301615-	15"	W25W301615-	16 ³ / ₄ "
16"	W25M301618-	18"	W25W301618-	16 ³ / ₄ "
16"	W25M301621-	21"	W25W301621-	16 ³ / ₄ "
16"	W25M301624-	24"	W25W301624-	16 ³ / ₄ "



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	W25M301330-	30"	W25W301230-	12 ³ / ₄ "
12 ³ / ₄ "	W25M301336-	36"	W25W301236-	12 ³ / ₄ "
12 ³ / ₄ "	W25M301342-	42"	W25W301242-	12 ³ / ₄ "
12 ³ / ₄ "	W25M301348-	48"	W25W301248-	12 ³ / ₄ "
16"	W25M301630-	30"	W25W301630-	16 ³ / ₄ "
16"	W25M301636-	36"	W25W301636-	16 ³ / ₄ "
16"	W25M301642-	42"	W25W301642-	16 ³ / ₄ "
16"	W25M301648-	48"	W25W301648-	16 ³ / ₄ "



Full Height Storage Cabinets – 84" High



Specifications:


Full Height Storage Cabinets are 84" high and 12", 16" or 22" deep. Lengths as shown.

Full Height Storage Cabinets are manufactured with four full-depth adjustable shelves and one fixed center shelf.



See page 95 for designating cabinet door style and door and cabinet material.

CAUTION: Full Height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall, Structural Module, or other cabinets to improve stability and ensure smooth door movement.


Open Storage Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S00M841318-0000	18"	S00W841218-00_0	12"
12 ³ / ₄ "	S00M841324-0000	24"	S00W841224-00_0	12"
16"	S00M841618-0000	18"	S00W841618-00_0	16"
16"	S00M841624-0000	24"	S00W841624-00_0	16"
22"	S00M842218-0000	18"	S00W842218-00_0	21 ³ / ₄ "
22"	S00M842224-0000	24"	S00W842224-00_0	21 ³ / ₄ "



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S00M841330-0000	30"	S00W841230-00_0	12"
12 ³ / ₄ "	S00M841336-0000	36"	S00W841236-00_0	12"
12 ³ / ₄ "	S00M841342-0000	42"	S00W841242-00_0	12"
12 ³ / ₄ "	S00M841348-0000	48"	S00W841248-00_0	12"
16"	S00M841630-0000	30"	S00W841630-00_0	16"
16"	S00M841636-0000	36"	S00W841636-00_0	16"
16"	S00M841642-0000	42"	S00W841642-00_0	16"
16"	S00M841648-0000	48"	S00W841648-00_0	16"
22"	S00M842230-0000	30"	S00W842230-00_0	21 ³ / ₄ "
22"	S00M842236-0000	36"	S00W842236-00_0	21 ³ / ₄ "
22"	S00M842242-0000	24"	S00W842242-00_0	21 ³ / ₄ "
22"	S00M842248-0000	24"	S00W842248-00_0	21 ³ / ₄ "



Full Height Storage Cabinets – 84" High



Specifications:

Full Height Storage Cabinets are 84" high and 12", 16" or 22" deep. Lengths as shown.


Full Height Storage Cabinets are manufactured with four full-depth adjustable shelves and one fixed center shelf.

Sliding Glazed and Panel Doors are suspended by rollers from a track at the top of the cabinet.


See page 95 for designating cabinet door style and door and cabinet material.

CAUTION: Full Height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall, Structural Module, or other cabinets to improve stability and ensure smooth door movement.


Glazed Sliding Door Cabinets




Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S11M841330-___0	30"	S11W841230-00_0	12"
12 ³ / ₄ "	S11M841336-___0	36"	S11W841236-00_0	12"
12 ³ / ₄ "	S11M841342-___0	42"	S11W841242-00_0	12"
12 ³ / ₄ "	S11M841348-___0	48"	S11W841248-00_0	12"
16"	S11M841630-___0	30"	S11W841630-00_0	16"
16"	S11M841636-___0	36"	S11W841636-00_0	16"
16"	S11M841642-___0	42"	S11W841642-00_0	16"
16"	S11M841648-___0	48"	S11W841648-00_0	16"
22"	S11M842230-___0	30"	S11W842230-00_0	21 ³ / ₄ "
22"	S11M842236-___0	36"	S11W842236-00_0	21 ³ / ₄ "
22"	S11M842242-___0	24"	S11W842242-00_0	21 ³ / ₄ "
22"	S11M842248-___0	24"	S11W842248-00_0	21 ³ / ₄ "




Paneled Sliding Door Cabinets



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S12M841330-___0	30"	S12W841230-00_0	12"
12 ³ / ₄ "	S12M841336-___0	36"	S12W841236-00_0	12"
12 ³ / ₄ "	S12M841342-___0	42"	S12W841242-00_0	12"
12 ³ / ₄ "	S12M841348-___0	48"	S12W841248-00_0	12"
16"	S12M841630-___0	30"	S12W841630-00_0	16"
16"	S12M841636-___0	36"	S12W841636-00_0	16"
16"	S12M841642-___0	42"	S12W841642-00_0	16"
16"	S12M841648-___0	48"	S12W841648-00_0	16"
22"	S12M842230-___0	30"	S12W842230-00_0	21 ³ / ₄ "
22"	S12M842236-___0	36"	S12W842236-00_0	21 ³ / ₄ "
22"	S12M842242-___0	24"	S12W842242-00_0	21 ³ / ₄ "
22"	S12M842248-___0	24"	S12W842248-00_0	21 ³ / ₄ "



Full Height Storage Cabinets – 84" High



Specifications:


Full Height Storage Cabinets are 84" high and 12", 16" or 22" deep. Lengths as shown.

Full Height Storage Cabinets are manufactured with four full-depth adjustable shelves and one fixed center shelf.

See page 95 for designating cabinet door style, door and cabinet material, and pull and hinge selection.

CAUTION: Full Height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall, Structural Module, or other cabinets to improve stability and ensure smooth door movement.

Glazed Swinging Door Cabinets





Left-Hand Swing Door


Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S20M841318L	18"	S20W841218L	12 ³ / ₄ "
12 ³ / ₄ "	S20M841324L	24"	S20W841224L	12 ³ / ₄ "
16"	S20M841618L	18"	S20W841618L	16 ³ / ₄ "
16"	S20M841624L	24"	S20W841624L	16 ³ / ₄ "
22"	S20M842218L	18"	S20W842218L	22"
22"	S20M842224L	24"	S20W842224L	22"

Right-Hand Swing Door


Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S20M841318-	18"	S20W841218-	12 ³ / ₄ "
12 ³ / ₄ "	S20M841324-	24"	S20W841224-	12 ³ / ₄ "
16"	S20M841618-	18"	S20W841618-	16 ³ / ₄ "
16"	S20M841624-	24"	S20W841624-	16 ³ / ₄ "
22"	S20M842218-	18"	S20W842218-	22"
22"	S20M842224-	24"	S20W842224-	22"



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S20M841330-	30"	S20W841230-	12 ³ / ₄ "
12 ³ / ₄ "	S20M841336-	36"	S20W841236-	12 ³ / ₄ "
12 ³ / ₄ "	S20M841342-	42"	S20W841242-	12 ³ / ₄ "
12 ³ / ₄ "	S20M841348-	48"	S20W841248-	12 ³ / ₄ "
16"	S20M841630-	30"	S20W841630-	16 ³ / ₄ "
16"	S20M841636-	36"	S20W841636-	16 ³ / ₄ "
16"	S20M841642-	42"	S20W841642-	16 ³ / ₄ "
16"	S20M841648-	48"	S20W841648-	16 ³ / ₄ "
22"	S20M842230-	30"	S20W842230-	21 ³ / ₄ "
22"	S20M842236-	36"	S20W842236-	21 ³ / ₄ "
22"	S20M842242-	24"	S20W842242-	21 ³ / ₄ "
22"	S20M842248-	24"	S20W842248-	21 ³ / ₄ "



Full Height Storage Cabinets – 84" High



Specifications:


Full Height Storage Cabinets are 84" high and 12", 16" or 22" deep. Lengths as shown.

Full Height Storage Cabinets are manufactured with four full-depth adjustable shelves and one fixed center shelf.

See page 95 for designating cabinet door style, door and cabinet material, and pull and hinge selection.

CAUTION: Full Height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall, Structural Module, or other cabinets to improve stability and ensure smooth door movement.

Paneled Swinging Door Cabinets





Left-Hand Swing Door


Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S25M841318L	18"	S25W841218L	12 ³ / ₄ "
12 ³ / ₄ "	S25M841324L	24"	S25W841224L	12 ³ / ₄ "
16"	S25M841618L	18"	S25W841618L	16 ³ / ₄ "
16"	S25M841624L	24"	S25W841624L	16 ³ / ₄ "
22"	S25M842218L	18"	S25W842218L	22"
22"	S25M842224L	24"	S25W842224L	22"

Right-Hand Swing Door


Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S25M841318-	18"	S25W841218-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841324-	24"	S25W841224-	12 ³ / ₄ "
16"	S25M841618-	18"	S25W841618-	16 ³ / ₄ "
16"	S25M841624-	24"	S25W841624-	16 ³ / ₄ "
22"	S25M842218-	18"	S25W842218-	22"
22"	S25M842224-	24"	S25W842224-	22"



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S25M841330-	30"	S25W841230-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841336-	36"	S25W841236-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841342-	42"	S25W841242-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841348-	48"	S25W841248-	12 ³ / ₄ "
16"	S25M841630-	30"	S25W841630-	16 ³ / ₄ "
16"	S25M841636-	36"	S25W841636-	16 ³ / ₄ "
16"	S25M841642-	42"	S25W841642-	16 ³ / ₄ "
16"	S25M841648-	48"	S25W841648-	16 ³ / ₄ "
22"	S25M842230-	30"	S25W842230-	21 ³ / ₄ "
22"	S25M842236-	36"	S25W842236-	21 ³ / ₄ "
22"	S25M842242-	24"	S25W842242-	21 ³ / ₄ "
22"	S25M842248-	24"	S25W842248-	21 ³ / ₄ "



Full Height Storage Cabinets – 84" High



Specifications:


Full Height Storage Cabinets are 84" high and 12", 16" or 22" deep. Lengths as shown.

Full Height Storage Cabinets are manufactured with four full-depth adjustable shelves and one fixed center shelf.

See page 95 for designating cabinet door style, door and cabinet material, and pull and hinge selection.

CAUTION: Full Height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall, Structural Module, or other cabinets to improve stability and ensure smooth door movement.

Paneled Swinging Door Cabinets





Left-Hand Swing Door

Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S25M841318L	18"	S25W841218L	12 ³ / ₄ "
12 ³ / ₄ "	S25M841324L	24"	S25W841224L	12 ³ / ₄ "
16"	S25M841618L	18"	S25W841618L	16 ³ / ₄ "
16"	S25M841624L	24"	S25W841624L	16 ³ / ₄ "
22"	S25M842218L	18"	S25W842218L	22"
22"	S25M842224L	24"	S25W842224L	22"

Right-Hand Swing Door

Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S25M841318-	18"	S25W841218-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841324-	24"	S25W841224-	12 ³ / ₄ "
16"	S25M841618-	18"	S25W841618-	16 ³ / ₄ "
16"	S25M841624-	24"	S25W841624-	16 ³ / ₄ "
22"	S25M842218-	18"	S25W842218-	22"
22"	S25M842224-	24"	S25W842224-	22"



Depth	Steel Cabinet	Length	Wood Cabinet	Depth
12 ³ / ₄ "	S25M841330-	30"	S25W841230-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841336-	36"	S25W841236-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841342-	42"	S25W841242-	12 ³ / ₄ "
12 ³ / ₄ "	S25M841348-	48"	S25W841248-	12 ³ / ₄ "
16"	S25M841630-	30"	S25W841630-	16 ³ / ₄ "
16"	S25M841636-	36"	S25W841636-	16 ³ / ₄ "
16"	S25M841642-	42"	S25W841642-	16 ³ / ₄ "
16"	S25M841648-	48"	S25W841648-	16 ³ / ₄ "
22"	S25M842230-	30"	S25W842230-	21 ³ / ₄ "
22"	S25M842236-	36"	S25W842236-	21 ³ / ₄ "
22"	S25M842242-	24"	S25W842242-	21 ³ / ₄ "
22"	S25M842248-	24"	S25W842248-	21 ³ / ₄ "

