

# Enterprise



**ADAPTABLE  
MOVEABLE WORKSTATIONS**

**KEWAUNEE<sup>®</sup>**

*...encouraging new discovery*

# ENTERPRISE

## *Freestanding Tables*

- Self-supporting workstations
- Single and double sided
- Adjustable height worksurfaces
- Integrated shelving systems
- Pre-plumbed integrated services
- Pre-wired integrated power and data





## Contents

Enterprise Workstations.....	2
Single-sided Support Frames.....	4
Double-sided Support Frames....	8
Worksurface Frames .....	12
Free-standing Tables .....	13
Work Tops .....	14
Enterprise Shelving.....	16
Steel & Wood Cabinets.....	22
Wood Mobile Cabinets.....	26
Steel Mobile Cabinets .....	28
Add-A-Drawer Cabinets.....	30
Mobile Cabinet Tops .....	31
Suspended Base Cabinets .....	32
Suspended Upper Cabinets .....	36
Ceiling Service Panels .....	38
Technical Data .....	42
Index.....	50



Kewaunee's Enterprise is a self-supporting moveable workstation table design. Enterprise is available with an above counter support structure, or as an adjustable height four legged table. Available pre-plumbed and pre-wired, this design easily incorporates with overhead carrier and panel systems. Worksurfaces are adjustable from sitting height to standing height, and overhead shelving is adjustable on 1" increments. Enterprise is designed for mobile, suspended, or floor mounted casework.

Kewaunee's Enterprise provides a completely moveable laboratory casework system to meet the needs of a changing lab environment.